

mackworth

TX Advanced Ceiling Track Lift



Set Up Guide


Contents

1.0	Introduction	3
2.0	Safety Precautions	4
2.1	Required Tools and Equipment for Lift Set Up and Decommissioning.....	5
3.0	Lift Unboxing.....	6
4.0	Lift Installation	7
4.1	Ceiling Track Access	7
4.2	Installing the Charger	8
4.2.1	Installing Charging Dock Type 1 (Track Type 1)	8
4.2.2	Installing Charging Dock Type 2 (Track Type 2)	9
4.2.3	Installing Charging Dock Type 3 (Track Type 3, 4 and 5).....	10
4.2.4	Installing Charging Dock Type 4 (Track Type 1) – QRT Lift.....	11
4.2.5	Installing Charging Dock Type 5 (Track Type 3, 4 and 5) – QRT Lift.....	12
4.3	Installing the Lift	13
4.3.1	Installing TX Advanced - Powered Traverse & Return to Charge (RTC)	13
4.3.2	Installing TX Advanced – Constant Charge (Powered Traverse and Manual Traverse).....	14
4.3.3	Installing TX Advanced – Powered Auxiliary (Powered Traverse and Manual Traverse)	15
4.3.4	Installing TX Advanced – Quick Release Track (QRT)	16
5.0	Component Attachment	17
5.1	Attaching the Handset	17
5.2	Attaching the Carry Bar.....	18
6.0	Lift Start Up and Function Checks.....	19
6.1	Turning the Lift ON and OFF	19
6.2	Operating Controls.....	19
6.3	Traversing the Lift along the Track	20
6.4	Charging the Lift.....	20
7.0	Accessories.....	20
8.0	Decommissioning.....	21
8.1	Lift Removal	21
8.2	Lift Dismantling.....	21
8.3	Lift Disposal.....	21

1.0 Introduction

Welcome to the Mackworth TX Advanced Ceiling Track Lift Commissioning and Decommissioning Manual, this manual includes the TX440 Advanced and the TX600 Advanced variants. Here it will explain on how to successfully install the lift into a pre-existing Ceiling Track safely and securely, as well as how to remove the lift from the Track. Use the contents page for easier navigation around the manual to discover the required section. The aim of the manual is that to help and assist you, so that you can provide a successful installation and removal of the lift for the satisfaction and safety of the customer.

The TX Advanced is available in different configurations dependent on track profile. The track profiles are shown below:



When Commissioning the Mackworth TX Advanced Lift, one of the above track profiles should match the track the lift is about to be installed into.

To ensure the lift available is compatible with the pre-installed track, refer to the table below.

This following guide covers the following TX Advanced Models:

TX400 Lift Type	Track Type	TX600 Lift Type	Track Type
TX440 Advanced PT – 122617	Type 1	TX600 Advanced PT – 122622	Type 1
TX440 Advanced PT – 122656	Type 3	TX600 Advanced PT – 122666	Type 3
TX440 Advanced PT – 122659	Type 4	TX600 Advanced PT – 122669	Type 4
TX440 Advanced PT – 122662	Type 5	TX600 Advanced PT – 122672	Type 5
TX440 Advanced QRT – 122678	Type 1	TX600 Advanced QRT – 122679	Type 1
TX440 Advanced QRT – 122688	Type 3	TX600 Advanced QRT – 122689	Type 3
TX440 Advanced QRT – 122693	Type 4	TX600 Advanced QRT – 122694	Type 4
TX440 Advanced QRT – 122698	Type 5	TX600 Advanced QRT – 122699	Type 5
TX440 Advanced RTC – 122703	Type 1	TX600 Advanced RTC – 122704	Type 1
TX440 Advanced RTC – 122708	Type 3	TX600 Advanced RTC – 122709	Type 3
TX440 Advanced RTC – 122713	Type 4	TX600 Advanced RTC – 122714	Type 4
TX440 Advanced RTC – 122718	Type 5	TX600 Advanced RTC – 122719	Type 5
TX440 Advanced MTCC – 122616	Type 1	TX600 Advanced MTCC – 122621	Type 1
TX440 Advanced MTCC – 122655	Type 2	TX600 Advanced MTCC – 122665	Type 2
TX440 Advanced MTCC – 122658	Type 3	TX600 Advanced MTCC – 122668	Type 3
TX440 Advanced MTCC – 122661	Type 4	TX600 Advanced MTCC – 122671	Type 4
TX440 Advanced MTCC – 122664	Type 5	TX600 Advanced MTCC – 122674	Type 5
TX440 Advanced PTCC – 122619	Type 1	TX600 Advanced PTCC – 122624	Type 1
TX440 Advanced PTCC – 122657	Type 3	TX600 Advanced PTCC – 122667	Type 3
TX440 Advanced PTCC – 122660	Type 4	TX600 Advanced PTCC – 122670	Type 4
TX440 Advanced PTCC – 122663	Type 5	TX600 Advanced PTCC – 122673	Type 5
TX440 Advanced MTPA – 122615	Type 1	TX600 Advanced MTPA – 122620	Type 1
TX440 Advanced PTPA – 122618	Type 1	TX440 Advanced PTPA – 122623	Type 1

Table 1

2.0 Safety Precautions

Read and understand this manual in its entirety before installing the Ceiling Track Lift.

The information in this manual is important for the safety of anyone near the TX lift and must be read and understood to help prevent injuries.

The TX Lift has been specifically designed to be installed in both the professional and home health care environments.

No matter the environment, health and safety factors should be considered to ensure the safety and essential performance of the lift and to avoid unnecessary damage or injuries to people within the area of the lift.

Typical examples include radiated heat (e.g. from a heater or fire place), excessive moisture impacting electrical performance (e.g. from a bathroom or kitchen area) and the correct storage of the lift after use (e.g. handset position on the carry bar).

The lift is not intended to be used in environments where there are rapid changes in the environmental temperature and humidity during intended use.

Any commissioning or decommissioning of the lift that has not been carried out by an authorized Mackworth personnel will void the product warranty.

- Mackworth authorised personnel must be the only Engineers who Commission and Decommission the TX lift.
- The lift must be turned OFF during installation and removal for the safety of the Engineer and Customer.
- Always ensure suitable clearance when installing and removing a lift into/from a Ceiling Track.
- Ensure that all the below procedures are followed correctly as instructed in this manual.
- All listed tools and equipment stated in this manual must be used to safely commission and decommission the lift.
- The TX lift has been specifically designed to be installed in both the professional and home health care environments only and should not be installed elsewhere.



Unauthorised modifications on any Mackworth product may affect its safety. The manufacturer will not be held responsible for any accident, incident or deficiencies of performance that occur as a result of any unauthorised modification to its products.

2.1 Required Tools and Equipment for Lift Set Up and Decommissioning

When carrying out work on the Ceiling Track Lift you will require the following:

Tools Required

2mm T Bar Allen wrench
2.5mm T Bar Allen wrench
3mm T Bar Allen wrench
4mm T Bar Allen wrench
5mm T Bar Allen wrench
5.5mm combination spanner
8mm combination spanner
13mm combination spanner x2
No. 2 Pozi drive screwdriver
Small slotted screwdriver
External Circlip pliers 3-10mm
Long nose pliers
Side snips
Power Drill
Drill Piece set up to 8.5mm

Equipment Required

Step Ladder
Marker Pen
Cloth
Lint Free Cloth
The correct documentation
Service Stand

3.0 Lift Unboxing

This section will summarize the layout of the Lift Packaging and what is included in the Box. It is recommended a small pocketknife is used for smoother unpacking of the lift. The lift is packed into a single box (11"x26.5"x14.5"), weighing approximately 24 lb (11 kg). Below is a list of components that are included in the box.

Item	Description	Product Code
1	TX Lift	See Table 1:
2	Carry Bar	300024
3	Handset	122073, 122081, 122349, 122075
4	Lift Charger	122100, 122326, 122327
5	Allen Key – 4mm	122093
6	User Manual	999553

Table 2




Figure 1

Using a knife to open the box down the middle, the box should open to a similar view to the image on the right. It will include all the above components.

Below is an image of all the included components for easier reference.

Remove the components from the box, taking the parts out in the following order:



Figure 2

1. Carry bar
2. Packaging foam
3. Handset
4. Lift battery charger
5. Trolley
6. User Manual
7. TX Lift

Once removed from the box, remove each individual component from their plastic bag ready for Track Installation.




Figure 3

4.0 Lift Installation

This section will cover the following guidelines on how to remove the safety components from the pre-installed track, how to install the charger into the track and how to install the lift into the track.

This section will cover all variant types of the TX lift range.

4.1 Ceiling Track Access

Prior to installing the Trolley into the pre-installed track, the track entrance must first be opened. A pre-installed track should include an end stop, safety bolt and an end cap. If these parts are not installed, it is essential that they are once the Trolley is installed.

Follow the guidance below to remove the ceiling track components. It is recommended that a step ladder is used to perform the task safer and easier.

1. Remove the end cap from the end of the track. (Figure 4)
2. Using two 13mm spanners, loosen the safety bolt and remove the securing nut from the end. (Figure 5)
3. Remove the safety bolt from the track and remove the large rubber bumper held within the track. (Figure 6)
4. Using a 5mm Allen Key, release the two screws securing the end stop and slide the end stop out of the track. (Figure 7)

To refit the end stop, safety bolt and end cap, follow the steps above in reverse order noting the following points:

- Ensure the end stop bumper is facing into the ceiling track.
- Ensure all fasteners are secure and cannot come undone without the use of tools.
- When fitting a charging dock, an end stop is not required at that end of the track.

(The images in this section all refer to a Track Type 1 installation, but the same process applies for all the remaining track types)



Figure 4



Figure 5



Figure 6



Figure 7

4.2 Installing the Charger

The section below will cover the Installation of the lift chargers into their specified Tracks. The list below specifies which chargers are compatible with which track types.

- Standard Charging Dock (122100) is designed for the Standard Prism Track (Track Type 1).
- Charging Dock Type 2 (122327) is designed for Track Type 2.
- Charging Dock Type 3 (122326) is designed for Track Type 3, 4 and 5.
- Charging Dock Type 4 (122338) is designed for Track Type 1 (QRT Lift)
- Charging Dock Type 5 (122339) is designed for Track Type 3, 4 and 5. (QRT Lift)

Use the Track Profile diagrams and the in the introduction to ensure that the lift and Charger provided is correct and compatible with the Pre-Installed Track. Do not use a variant lift and charger with an incorrect track as they are incompatible.

4.2.1 Installing Charging Dock Type 1 (Track Type 1)

To install the TX Charger, follow the guidance below:

1. Complete all relevant steps in 4.1 to gain access to the Ceiling Track.
2. A step ladder is recommended for easier and safer installation.
3. Align the Charger as shown in Figure 8, the Charger must be inserted into the track backward (as shown) and must slid all the way to the other end of the Track. This ensures that the Charger is facing inward and allow the lift to dock.
4. Position the Charger in front of the Safety Bolt bumper to achieve maximum track length. (Figure 9)
5. Fix the Charger in place using a 5mm Allen Key to tighten the two bolts. Ensure the Charger is positioned straight in the track. (Figure 10)
6. Route the Charger wire in a tidy fashion to the closest power supply.
7. The Charger block locates into its bracket as shown in figure 11. Raw Plugs and Screws are provided to fix the unit to the wall if desired.
8. Choose a suitable location on the buildings wall where the charging unit can be fitted. The wall must be of good structure to fit the three screw fixing points.
9. Drill three 7mm holes into the wall that will allow the three given brown raw plugs to be fitted, measure the gap between the three drilled holes to match the profile of the of the charging unit bracket.
10. Place the raw plugs into the drilled holes.
11. Align the charging bracket fixing point with the holes and use a pozi head screwdriver to fix the two screws and secure the bracket to the wall. (Figure 12)
12. Fit the Charging Unit into the Bracket
13. Plug the Charger into the nearest wall socket and turn on the power supply
14. The Charging Dock is ready for use.
15. If a lift is not being installed directly after charger installation. Refer back to section 4.1 to refit the removed track components to close access to the Ceiling Track.



Figure 8



Figure 9



Figure 10



Figure 11



Figure 12

4.2.2 Installing Charging Dock Type 2 (Track Type 2)

Although the charger looks different to the Type 1 Charging Dock, the installation process is the same. See below for details on installing the variant along with the variant charger images.

To install the TX Charger, follow the guidance below:

1. Complete all relevant steps in 4.1 to gain access to the Ceiling Track.
2. A step ladder is recommended for easier and safer installation.
3. Align the Charger as shown in Figure 13, the Charger must be inserted into the track backward (as shown) and must slid all the way to the other end of the Track. This ensures that the Charger is facing inward and allow the lift to dock.
4. Position the Charger in front of the Safety Bolt bumper to achieve maximum track length. (Figure 14)
5. Fix the Charger in place using a 5mm Allen Key to tighten the two bolts. Ensure the Charger is positioned straight in the track. (Figure 15)
6. Route the Charger wire in a tidy fashion to the closest power supply.
7. The Charger block locates into its bracket as shown in figure 16. Raw Plugs and Screws are provided to fix the unit to the wall if desired.
8. Choose a suitable location on the buildings wall where the charging unit can be fitted. The wall must be of good structure to fit the three screw fixing points.
9. Drill three 7mm holes into the wall that will allow the three given brown raw plugs to be fitted, measure the gap between the three drilled holes to match the profile of the of the charging unit bracket.
10. Place the raw plugs into the drilled holes.
11. Align the charging bracket fixing point with the holes and use a pozi head screwdriver to fix the two screws and secure the bracket to the wall. (Figure 17)
12. Fit the Charging Unit into the Bracket
13. Plug the Charger into the nearest wall socket and turn on the power supply
14. The Charging Dock is ready for use.
15. If a lift is not being installed directly after charger installation. Refer back to section 4.1 to refit the removed track components to close access to the Ceiling Track.



Figure 13



Figure 14



Figure 15



Figure 16



Figure 17

4.2.3 Installing Charging Dock Type 3 (Track Type 3, 4 and 5)

Although the charger looks different to the Type 1 Charging Dock, the installation process is the same. See below for details on installing the variant along with the variant charger images.

To install the TX Charger, follow the guidance below:

1. Complete all relevant steps in 4.1 to gain access to the Ceiling Track.
2. A step ladder is recommended for easier and safer installation.
3. Align the Charger as shown in Figure 18, the Charger must be inserted into the track backward (as shown) and must slid all the way to the other end of the Track. This ensures that the Charger is facing inward and allow the lift to dock.
4. Position the Charger in front of the Safety Bolt bumper to achieve maximum track length. (Figure 19)
5. Fix the Charger in place using a 5mm Allen Key to tighten the two bolts. Ensure the Charger is positioned straight in the track. (Figure 20)
6. Route the Charger wire in a tidy fashion to the closest power supply.
7. The Charger block locates into its bracket as shown in figure 21. Raw Plugs and Screws are provided to fix the unit to the wall if desired.
8. Choose a suitable location on the buildings wall where the charging unit can be fitted. The wall must be of good structure to fit the three screw fixing points.
9. Drill three 7mm holes into the wall that will allow the three given brown raw plugs to be fitted, measure the gap between the three drilled holes to match the profile of the of the charging unit bracket.
10. Place the raw plugs into the drilled holes.
11. Align the charging bracket fixing point with the holes and use a pozi head screwdriver to fix the two screws and secure the bracket to the wall. (Figure 22)
12. Fit the Charging Unit into the Bracket
13. Plug the Charger into the nearest wall socket and turn on the power supply
14. The Charging Dock is ready for use.
15. If a lift is not being installed directly after charger installation. Refer back to section 4.1 to refit the removed track components to close access to the Ceiling Track.




Figure 18



Figure 19



Figure 20



Figure 21



Figure 22

4.2.4 Installing Charging Dock Type 4 (Track Type 1) - QRT Lift

To install the TX Charger, follow the guidance below:

1. Complete all relevant steps in 4.1 to gain access to the Ceiling Track.
2. A step ladder is recommended for easier and safer installation.
3. Align the Charger as shown in Figure 23, the Charger must be inserted into the track backward (as shown) and must slid all the way to the other end of the Track. This ensures that the Charger is facing inward and allow the lift to dock.
4. Position the Charger in front of the Safety Bolt bumper to achieve maximum track length. (Figure 24)
5. Fix the Charger in place using a 5mm Allen Key to tighten the two bolts. Ensure the Charger is positioned straight in the track. (Figure 25)
6. Place a self-adhesive clip onto the track directly above the Charger. (Figure 26)
7. Fit the locking pin attached to the charger via a chain into the clip for safe keeping (Figure 26)
8. Route the Charger wire in a tidy fashion to the closest power supply.
9. The Charger block locates into its bracket as shown in figure 27. Raw Plugs and Screws are provided to fix the unit to the wall if desired.
10. Choose a suitable location on the buildings wall where the charging unit can be fitted. The wall must be of good structure to fit the three screw fixing points.
11. Drill three 7mm holes into the wall that will allow the three given brown raw plugs to be fitted, measure the gap between the three drilled holes to match the profile of the of the charging unit bracket.
12. Place the raw plugs into the drilled holes.
13. Align the charging bracket fixing point with the holes and use a pozi head screwdriver to fix the two screws and secure the bracket to the wall. (Figure 28)
14. Fit the Charging Unit into the Bracket.
15. Plug the Charger into the nearest wall socket and turn on the power supply.
16. Drill a 9mm hole in front of the charger for the trolley set dismount point.
17. Place the trolley set dismount point sticker onto the hole as shown in figure 29.
18. The Charging Dock is ready for use.
19. If a lift is not being installed directly after charger installation. Refer back to section 4.1 to refit the removed track components to close access to the Ceiling Track.




Figure 23



Figure 24



Figure 25



Figure 26



Figure 27



Figure 28



Figure 29

4.2.5 Installing Charging Dock Type 5 (Track Type 3, 4 and 5) - QRT Lift

To install the TX Charger, follow the guidance below:

1. Complete all relevant steps in 4.1 to gain access to the Ceiling Track.
2. A step ladder is recommended for easier and safer installation.
3. Align the Charger as shown in Figure 30, the Charger must be inserted into the track backward (as shown) and must slid all the way to the other end of the Track. This ensures that the Charger is facing inward and allow the lift to dock.
4. Position the Charger in front of the Safety Bolt bumper to achieve maximum track length. (Figure 31)
5. Fix the Charger in place using a 5mm Allen Key to tighten the two bolts. Ensure the Charger is positioned straight in the track. (Figure 32)
6. Place a self-adhesive clip onto the track directly above the Charger. (Figure 33)
7. Fit the locking pin attached to the charger via a chain into the clip for safe keeping (Figure 33)
8. Route the Charger wire in a tidy fashion to the closest power supply.
9. The Charger block locates into its bracket as shown in figure 34. Raw Plugs and Screws are provided to fix the unit to the wall if desired.
10. Choose a suitable location on the buildings wall where the charging unit can be fitted. The wall must be of good structure to fit the three screw fixing points.
11. Drill three 7mm holes into the wall that will allow the three given brown raw plugs to be fitted, measure the gap between the three drilled holes to match the profile of the of the charging unit bracket.
12. Place the raw plugs into the drilled holes.
13. Align the charging bracket fixing point with the holes and use a pozi head screwdriver to fix the two screws and secure the bracket to the wall. (Figure 35)
14. Fit the Charging Unit into the Bracket.
15. Plug the Charger into the nearest wall socket and turn on the power supply.
16. Drill a 9mm hole in front of the charger for the trolley set dismount point.
17. Place the trolley set dismount point sticker onto the hole as shown in figure 36.
18. The Charging Dock is ready for use.
19. If a lift is not being installed directly after charger installation. Refer back to section 4.1 to refit the removed track components to close access to the Ceiling Track.




Figure 30



Figure 31



Figure 32



Figure 33



Figure 34



Figure 35




Figure 36

4.3 Installing the Lift

The section below will cover the Installation of the lifts into their specified Tracks. Refer to the Table in the introduction to determine if the lift available is compatible with the pre-installed track. Do not use a variant lift and track together as they are incompatible.

4.3.1 Installing TX Advanced - Powered Traverse & Return to Charge (RTC)

To install the TX lift, follow the guidance below:

1. Complete all the relevant steps in section 4.1.
2. Complete all the relevant steps in section 4.2.
3. A step ladder is recommended for easier and safer installation.
4. Using both hands to support the lift, align the lift horizontally with the Ceiling Track entrance. (Figure 37)
5. Ensure the charging port is **leading** into the track. The charging port is circled in (Figure 37)
6. Carefully place the lift wheel assembly into the track without damaging the lift components. (Figure 38)
7. Using the wheels, slide the lift away from the end of the ceiling track for safety. (Figure 39 and 40)
8. Reverse the Steps in section 4.1 to secure the ceiling track entrance. Fitting the end stop, safety bolt and end cap correctly.
9. Lift Installation Complete.

To remove the lift, follow the steps above in reverse.



Figure 37



Figure 38



Figure 39



Figure 40

4.3.2 Installing TX Advanced - Constant Charge (Powered Traverse and Manual Traverse)

To install the TX lift, follow the guidance below:

1. Complete all the relevant steps in section 4.1.
2. Complete all the relevant steps in section 4.2.
3. A step ladder is recommended for easier and safer installation.
4. Using both hands to support the lift, align the lift horizontally with the Ceiling Track entrance.
5. Ensure the positive symbol on the constant charge port is facing the same direction as the positive symbol on the pre-installed track. The symbols are circled in (Figure 41)
6. Carefully place the lift wheel assembly into the track without damaging the lift components. (Figure 42 and 43)
7. Using the wheels, slide the lift away from the end of the ceiling track for safety. (Figure 44 and 45)
8. Reverse the Steps in section 4.1 to secure the ceiling track entrance. Fitting the end stop, safety bolt and end cap correctly.
9. Left Installation Complete.

To remove the lift, follow the steps above in reverse.




Figure 41



Figure 42




Figure 43



Figure 44



Figure 45

4.3.3 Installing TX Advanced - Powered Auxiliary (Powered Traverse and Manual Traverse)

To install the TX lift, follow the guidance below:

1. Complete all the relevant steps in section 4.1.
2. Complete all the relevant steps in section 4.2.
3. A step ladder is recommended for easier and safer installation.
4. Using both hands to support the lift, align the lift horizontally with the Ceiling Track entrance. (Figure 46)
5. Ensure the charging port is **leading** into the track. The charging port is circled in (Figure 47)
6. Carefully place the lift wheel assembly into the track without damaging the lift components. (Figure 48)
7. Using the wheels, slide the lift away from the end of the ceiling track for safety. (Figure 49)
8. Reverse the Steps in section 4.1 to secure the ceiling track entrance. Fitting the end stop, safety bolt and end cap correctly.
9. Left Installation Complete.

To remove the lift, follow the steps above in reverse.



Figure 46



Figure 47



Figure 48



Figure 49

4.3.4 Installing TX Advanced - Quick Release Track (QRT)

For a TX with a QRT system, the trolley must first be installed into the track before attaching the lift to the trolley. To install the TX Trolley, follow the guidance below:

1. Complete all the relevant steps in section 4.1.
2. Complete all the relevant steps in section 4.2.
3. A step ladder is recommended for easier and safer installation.
4. Align the trolley horizontally with the Ceiling Track entrance, with the attachment L-Shaped hooks facing downward. (Figure 50)
5. Ensure the charging port is **leading** into the track. The charging port is circled in (Figure 50)
6. Carefully place the Trolley into the track without damaging the component. (Figure 51)
7. Using the wheels, slide the Trolley away from the end of the ceiling track for safety. (Figure 52 and 53)
8. Reverse the Steps in section 4.1 to secure the ceiling track entrance. Fitting the end stop, safety bolt and end cap correctly.
9. Trolley Installation Complete.

To remove the Trolley, follow the steps above in reverse.



Figure 50



Figure 51



Figure 52



Figure 53

To install the lift onto the Trolley, follow the guidance below:

1. Complete all the relevant steps in section 4.1.
2. Complete all the relevant steps in 4.2.
3. Complete all the relevant steps in 4.3.4.
4. A step ladder is recommended for easier and safer installation.
5. Disengage the plunger on the trolley to allow attachment access.
6. Raise the TX to the height of the Trolley as shown. (Figure 54)
7. Align the anchor points with the trolley hooks, into its profile. (Figure 55)
8. Engage the plunger to lock the lift into the Trolley.
9. Plug the Charger lead into the lift cover. (Figure 56)
10. Cable tie the lead onto the chassis frame as shown. (Figure 57)
11. Lift installation complete.

To remove the lift, follow the steps above in reverse order noting the following points:

- Traverse the lift to the Trolley set dismount point, ensure to align the hole in the trolley with the hole in the track.
- Remove the Pin from the clip and place it through the track to secure the Trolley, this will allow the lift to be dismantled. Then follow the installation guide in reverse to detach the lift from the trolley.




Figure 54



Figure 55




Figure 56




Figure 57

5.0 Component Attachment

The section below will cover the attachment of the Carry Bar and Handset to the lift, as well as other essential requirements to fully commission the lift for use. The below procedure is relevant for both the TX600 Advanced and TX440 Advanced.



A sturdy ladder or steps may be required in order to access the underside of the lift to attach the handset.

5.1 Attaching the Handset

In this section the correct procedure on attaching the Handset will be determined.

The section refers to a four-control handset, but the same process applies to a two-way and six-way handset.

The Pneumatic Handset is attached to the lift through the Airline Grommet on the underside of the lift as shown in Figure 58. While holding the Airline Connector at the end of the Handset Air Tubes, locate the two small metal ribbed pins with the Airline Grommet on the lift. To correctly locate, ensure the profile of the two parts are lining up, this means the grey rib profile, which faces inwards on the lift aligns with the grey rib on the Handset, this is the only way the Handset can be connected, therefore it is not possible to incorrectly mate the two parts.

Once the two parts have been aligned, ensure that the Handset is pushed into the Grommet properly, with both metal pins being fully inserted into the Grommet. This is important to ensure there is no air leak which will affect the lift functionality.




Figure 58




Figure 59




Figure 60

To remove the handset, follow the procedure above in reverse.

5.2 Attaching the Carry Bar

Follow the procedure below to attach the Carry Bar to the Quick Release Hook located at the end of the Lift Tape.

1. Open the hooks locking mechanism by pressing the red tab downward as shown. (Figure 61 and 62)
2. Position the carry bar horizontally as shown in figure 63 and place the attachment point into the open hook.
3. Pivot the carry bar down to its natural position and release the hooks locking mechanism to secure the carry bar. (Figure 64)
4. Alternatively, with the carry bar positioned sideways, the pin at the top of the carry bar can gently push the locking mechanism out of the way as the pin is carefully moved into the hook, the locking mechanism will close itself.

To remove the Carry Bar, follow the steps above in reverse.



Figure 61



Figure 62




Figure 63




Figure 64

6.0 Lift Start Up and Function Checks

The section below will show how to start up the lift for the use of the customer as well as completing some basic functional checks to ensure the lift is in working order. The images below refer to the TX600 Advanced, but the same procedure applies to all models.

6.1 Turning the Lift ON and OFF

To start up the lift the Emergency Lower Tab must be located in the ON position, as the lift is delivered to the customer in the OFF position. To give power to the lift, the tab must be pressed upward vertically as shown in the image below. (Figure 65)

Press any button on the Handset once – This should start up the lift. The LCD Display should light up Green showing the Home Screen Display “Battery 100%”. (Figure 66). The LED should also light up green.



Figure 65





Figure 66

6.2 Operating Controls

By pressing the UP or the DOWN arrow button on the handset, the carry bar can be raised or lowered to the correct height for attaching the sling or positioning an individual. The UP/DOWN functions of the handset buttons are in relation to the travel of the lift. That is, the grey button at the top end of the handset activates the UP motion of the carry bar and the Green button activates the DOWN motion (Figure 67). These functions are found on all variants of the lift.

The Blue and Yellow buttons are used for Traversing the lift, this function is found on all powered traverse lifts. The colors will correspond with the arrows on the underside of the lift cover. Allowing to traverse one way or the other.

The Black and White buttons allow the lift to communicate with a Powered Auxiliary such as a turntable. These buttons are only available on the Powered Auxiliary option lifts. Activating these buttons when the lift is docked in a turntable will allow the turntable to rotate in the given direction.



- Up Function
- Down Function
- Traverse Left/Right
- Traverse Left/Right
- Powered Auxiliary Left/Right
- Powered Auxiliary Left/Right

Figure 67

6.3 Traversing the Lift along the Track

Always use extreme care when moving the lift along the track. Watch out for and avoid any obstructions that may cause injury to the individual in the sling, damage to the lift and/or to the obstruction. This applies to both manual and powered traversing.



After use, the lift should be located at the correct end of the track system for re-charging.

When needed, the lift should be moved along the track to a position directly above the person using the following appropriate method:

The operator shall lower the carry bar to a comfortable height so that it can easily be held by hand.

Gently pushing the carry bar, move the lift along the track to the required position, or traverse using the handset, where the person can be safely loaded into the sling.



NEVER pull the lift along the track using the handset or Emergency Lowering Cord as this could have a detrimental effect on the performance of the lift.

6.4 Charging the Lift

The TX lift is charged by the charging dock installed in the track. A charging dock should have been fitted during the installation of the lift.

To Traverse the TX600 or TX440, simply manually/powerd traverse the lift along the track to the charging dock at the far end of the track. Or for RTC versions use the RTC function to return the lift to the charging dock. Once the lift has made contact, it should begin charging immediately. The LED will begin to flash orange when charging, or if fully charged, will show a steady green. The LCD screen will display 'charging' or 'charged'.

At the end of each use of the lift, it is recommended that it be returned to the charging dock for placement and charging. This will ensure that the batteries are charged on a regular basis for peak performance and maximum life expectancy.




Figure 67

7.0 Accessories

The Remaining parts in the box include the User Manual and the Allen key which is used for Manual Emergency Lowering of the lift. These should be given to the customer for safe keeping. The Box should also be kept with the User and stored in a safe place. The Box will be required for decommissioning and transportation.

8.0 Decommissioning

The following section will instruct on how to correctly decommission the TX lift at the end of use. This section is relevant to both the TX600 Advanced and TX440 Advanced.



The product may be contaminated and has to be disinfected before decommissioning. See section 'Cleaning' in the User Manual for details of how to do this.

8.1 Lift Removal

The removal of the lift is essentially the Commissioning Steps in reversal. A summary of the correct order on how to remove the TX and its components has been listed below. For in detailed information on each of the below steps, see the relevant section within Commissioning for assistance.

1. Turn OFF the lift by pulling down the E-Lower Cord to its middle switch.
2. Detach the Carry Bar from the lift.
3. Detach the Handset from the lift
4. At the end of the Track, remove the End Cap, Safety Bolt and End Stop.
5. Carefully slide the lift out of the Track.
6. Turn OFF the Charger at the Wall Socket and remove the plug.
7. Remove the charging unit from the wall using a pozi head screwdriver.
8. Detach and carefully slide the Charger out of the Track.
9. Ensure the End Stop, Safety Bolt and End Cap are refitted to the Track once all the above steps are complete.
10. Lift Removal is complete.

8.2 Lift Dismantling

All serviceable components of the lift should be dismantled prior to the disposal of the lift. For guidance on how to dismantle all of the lift components, refer to the TX Advanced Service Manual and follow each sub-section within Section 9.0 - Servicing – Removal and Replacement.

8.3 Lift Disposal



When the lift has completed its life cycle and can no longer perform to its intended use safely the lift must be decommissioned by an approved Service Engineer. The following specifies the importance of correct disposal procedure including local laws and being environmentally friendly.

Please observe the local laws on recycling and respect the current laws for disposal within the community the device is being used within. If there is any uncertainty of the below guidelines, contact your local authorities to determine the proper method of disposal of potentially biohazardous parts and accessories.

The relevant components utilized in the manufacture of the device that can be recycled at the end of the device life are:

Fully recyclables:	Consideration when Recycling:
Chassis	Batteries
Plastic Covers	Wiring Looms – electronics
Metallic Internals – Hub etc.	PCB
Initial packaging of the device (cardboard)	Hand Control
Metallic fixing – Screws etc.	Motors
Plastic Mouldings – Wheels etc.	Lift Tape
Carry Bar	Charger

Ensure that this list is used as guidance and that the local laws in the given community overrule the suggested component disposal in the table above.

Dealer/service contact details:

Manufacturer contact details:

Prism Medical UK

Address: Unit 1 • Tir Llwyd Industrial Estate • St Asaph Avenue • Kinmel Bay • Conwy • LL18 5JZ

Telephone Number: 01924 840100

Disclaimer

While every effort has been made to ensure the accuracy of information contained in this installation manual, no liability can be accepted by Mackworth for any errors or omissions. Mackworth operates a policy of continuous improvement. Specifications and other data are subject to change without notice.



mackworth

54 West Industrial Drive • O'Fallon, MO 63366, USA